

# pcDNA6/TR

**A regulatory vector designed for use with the T-REx™ System**

**Catalog no. V1025-20**

**Version D**

012502

25-0274



[www.invitrogen.com](http://www.invitrogen.com)  
[tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)



# Table of Contents

Table of Contents .....	iii
General Information .....	iv
<b>Methods.....</b>	<b>1</b>
Using pcDNA6/TR.....	1
Creation of Stable Cell Lines .....	4
<b>Appendix .....</b>	<b>6</b>
Blasticidin .....	6
pcDNA6/TR Vector .....	7
Technical Service .....	9
Purchaser Notification.....	11
References .....	14

## General Information

---

<b>Contents</b>	20 µg of pcDNA6/TR, lyophilized in TE, pH 8.0
<b>Shipping/Storage</b>	Lyophilized plasmid is shipped at room temperature and should be stored at -20°C.

---

# Methods

## Using pcDNA6/TR

### Introduction

pcDNA6/TR is a 6.7 kb vector designed for use with the T-REx™ System (Catalog nos. K1020-01, K1020-02, K1030-01, and K1030-02). The vector expresses high levels of the tetracycline (Tet) repressor under the control of the human cytomegalovirus immediate-early (CMV) promoter. High-level stable and transient expression of the Tet repressor can be carried out in most mammalian cells. Tetracycline-regulated expression of a gene of interest may then be tested by transfecting the inducible expression plasmid into host cells expressing the Tet repressor. A map of pcDNA6/TR and a description of the features of the vector may be found in the **Appendix**, pages 7-8.

For more information about the T-REx™ System, please refer to the T-REx™ System manual. For information about T-REx™ inducible expression vectors, please refer to the manual for each specific vector. Manuals are available for downloading from our Web site ([www.invitrogen.com](http://www.invitrogen.com)) or by contacting Technical Service (see page 9). To order components of the T-REx™ System separately, please see **Accessory Products** below.

### A Note about the TetR Gene

The *TetR* gene used in pcDNA6/TR was originally isolated from the Tn10 transposon which confers resistance to tetracycline in *E. coli* and other enteric bacteria (Postle *et al.*, 1984). The *TetR* gene from Tn10 encodes a class B Tet repressor and is often referred to as *TetR(B)* in the literature (Hillen and Berens, 1994).

The *TetR* gene encodes a repressor protein of 207 amino acids with a calculated molecular weight of 23 kDa. For more information about the Tet repressor, its interaction with Tet operator sequences, and tetracycline regulation, please refer to the T-REx™ System Manual or to published reviews (Hillen and Berens, 1994; Hillen *et al.*, 1983).

### Accessory Products

Many of the reagents used in the T-REx™ System are available separately from Invitrogen. Ordering information is provided below.

Item	Amount	Purpose	Catalog no.
pcDNA4/TO	20 µg, lyophilized	Inducible expression vector	V1020-20
pcDNA4/TO/ <i>myc</i> -His A, B, C	20 µg each, lyophilized	Inducible expression vector	V1030-20
Blasticidin	50 mg, powder	Selection agent for regulatory plasmid	R210-01
Tetracycline	5 g, powder	Inducing agent	Q100-19
Zeocin™	1 g	Selection agent for	R250-01
	5 g	inducible expression vector	R250-05

*continued on next page*

# Using pcDNA6/TR, continued

---

## Maintenance of pcDNA6/TR

The pcDNA6/TR vector contains the ampicillin resistance gene and the blasticidin resistance gene, either of which allows selection of the plasmid in *E. coli*. To propagate and maintain the pcDNA6/TR vector, we recommend resuspending the vector in 20  $\mu$ l sterile water to prepare a 1  $\mu$ g/ $\mu$ l stock solution. Store the stock solution at -20°C.

Use this stock solution to transform a *recA*, *endA* *E. coli* strain like TOP10F' (Catalog no. C615-00), DH5 $\alpha$ F', JM109, INV $\alpha$ F', or equivalent. Select transformants on LB agar plates containing either 50 to 100  $\mu$ g/ml ampicillin or 100  $\mu$ g/ml blasticidin in Low Salt LB (see recipe below).

---

## Selection in *E. coli*

To facilitate selection of blasticidin-resistant *E. coli*, the salt concentration of the medium must remain low (< 90 mM) and the pH must be 7.0. Prepare LB broth and plates using the recipe below.

**Failure to lower the salt content of your LB medium will result in non-selection due to inhibition of the drug unless a higher concentration of blasticidin is used.**

---

## Low Salt LB Medium with Blasticidin

### Low Salt LB Medium:

10 g Tryptone  
5 g NaCl  
5 g Yeast Extract

1. Combine the dry reagents above and add deionized, distilled water to 950 ml. Adjust pH to 7.0 with 1 N NaOH. Bring the volume up to 1 liter. For plates, add 15 g/L agar before autoclaving.
  2. Autoclave on liquid cycle at 15 psi and 121°C for 20 minutes.
  3. Allow the medium to cool to at least 55°C before adding the blasticidin to 100  $\mu$ g/ml final concentration.
  4. Store plates at +4°C in the dark. Plates containing blasticidin are stable for up to 2 weeks.
- 

## Preparing a Glycerol Stock

Once you have identified the correct clone, purify the colony and make a glycerol stock for long-term storage. It is also a good idea to keep a DNA stock of your plasmid at -20°C.

- Streak the original colony out on an LB plate containing 50-100  $\mu$ g/ml ampicillin or 100  $\mu$ g/ml blasticidin in Low Salt LB. Incubate the plate at 37°C overnight.
  - Isolate a single colony and inoculate into 1-2 ml of LB containing 50-100  $\mu$ g/ml ampicillin or 100  $\mu$ g/ml blasticidin in Low Salt LB.
  - Grow the culture to mid-log phase ( $OD_{600} = 0.5-0.7$ ).
  - Mix 0.85 ml of culture with 0.15 ml of sterile glycerol and transfer to a cryovial.
  - Store at -80°C.
- 

*continued on next page*

## Using pcDNA6/TR, continued

---

### Plasmid Preparation

Plasmid DNA for transfection into eukaryotic cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells, and salt will interfere with lipid complexing, decreasing transfection efficiency. We recommend isolating plasmid DNA using the S.N.A.P.<sup>™</sup> MiniPrep Kit (10-15 µg DNA, Catalog no. K1900-01), the S.N.A.P.<sup>™</sup> MidiPrep Kit (10-200 µg DNA, Catalog no. K1910-01), or CsCl gradient centrifugation.

---

### Methods of Transfection

To transfect established mammalian cell lines (e.g. HeLa, COS-1) with pcDNA6/TR, please consult original references or the supplier of your cell line for the optimal method of transfection. We recommend that you follow exactly the protocol for your cell line. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler *et al.*, 1977), lipid-mediated (Felgner *et al.*, 1989; Felgner and Ringold, 1989) and electroporation (Chu *et al.*, 1987; Shigekawa and Dower, 1988). Invitrogen offers the Calcium Phosphate Transfection Kit and Lipofectamine<sup>™</sup> 2000 Reagent for mammalian cell transfection. For more information on transfection reagents available from Invitrogen, refer to our Web site ([www.invitrogen.com](http://www.invitrogen.com)) or contact Technical Service (see page 9).

Catalog no.	Description	Quantity
K2780-01	Calcium Phosphate Transfection Kit	75 reactions
11668-019	Lipofectamine <sup>™</sup> 2000 Reagent	1.5 ml

---

### Transient Transfection

You may use any of the methods above to transiently cotransfect pcDNA6/TR and your inducible expression construct into the mammalian host cell line. Because the amount of Tet repressor expressed in the cell will determine the level of transcriptional repression of the hybrid CMV/TetO<sub>2</sub> promoter in the inducible expression plasmid, we recommend that you increase the amount of pcDNA6/TR DNA transfected into your host cell line relative to inducible expression plasmid DNA. Increasing the ratio of pcDNA6/TR:inducible expression plasmid DNA from 1:1 to **at least 6:1** should ensure that a sufficient amount of Tet repressor is expressed to suitably repress basal transcription of your gene of interest. For more information about transfection and induction of expression using tetracycline, please refer to the T-REx<sup>™</sup> System manual.

---

# Creation of Stable Cell Lines

---

## Introduction

To generate a cell line that stably expresses the Tet repressor, you may transfect pcDNA6/TR into your mammalian host cell line and select with blasticidin. Cell lines expressing suitably high levels of the Tet repressor may then be used as hosts to stably or transiently express your gene of interest from the inducible expression vector. Before transfection, we recommend that you first test the sensitivity of your mammalian host cell to blasticidin as natural resistance varies among cell lines.

---

## Determination of Antibiotic Sensitivity

To generate a stable cell line expressing pcDNA6/TR, you need to determine the minimum concentration of blasticidin required to kill your untransfected host cell line. Typically, concentrations between 2 and 10 µg/ml blasticidin are sufficient to kill the untransfected host cell line. Test a range of concentrations (see below) to ensure that you determine the minimum concentration necessary for your cell line. For more information about blasticidin and instructions for use, please refer to the **Appendix**, page 6.

- Plate or split a confluent plate so the cells will be approximately 25% confluent. Prepare 6 plates of cells.
  - The next day, substitute culture medium with medium containing varying concentrations of blasticidin (e.g. 0, 1, 3, 5, 7.5, and 10 µg/ml).
  - Replenish the selective medium every 3-4 days. Cells sensitive to blasticidin will round up and detach from the plate. Dead cells will accumulate in the medium.
  - Count the number of viable cells at regular intervals to determine the appropriate concentration of blasticidin that prevents growth within 1-2 weeks after addition of blasticidin.
- 

## Possible Sites for Linearization

To obtain stable transfectants, you may choose to linearize the pcDNA6/TR plasmid before transfection. While linearizing your vector may not improve the efficiency of transfection, it increases the chances that the vector does not integrate in a way that disrupts either the *TetR* gene or other elements required for mammalian expression. The table below lists unique sites that may be used to linearize pcDNA6/TR prior to transfection. **Other restriction sites are possible.**

Enzyme	Restriction Site (bp)	Location	Supplier
<i>Bst</i> 1107 I	4470	Backbone	AGS*, Fermentas, Takara
<i>Sap</i> I	4733	Backbone	New England Biolabs
<i>Bsp</i> LU11 I	4849	Backbone	Boehringer-Mannheim
<i>Eam</i> 1105 I	5739	Ampicillin gene	AGS*, Fermentas, Takara
<i>Fsp</i> I	5961	Ampicillin gene	Many

\*Angewandte Gentechnologie Systeme

---

*continued on next page*



# Creation of Stable Cell Lines, continued

---

## Selection of Stable Integrants

Once you have determined the appropriate blasticidin concentration to use for selection, you can generate a stable cell line expressing pcDNA6/TR.

1. Transfect your cell line of choice with pcDNA6/TR using the desired protocol. Include a sample of untransfected cells as a negative control.
  2. 24 hours after transfection, wash the cells and add fresh medium to the cells.
  3. 48 hours after transfection, split the cells into fresh medium containing blasticidin at the appropriate concentration for your cell line. Split the cells such that they are no more than 25% confluent. If the cells are too dense, the blasticidin will not kill the untransfected cells.
  4. Replenish selective medium every 3-4 days until blasticidin-resistant colonies are detected. Typically, blasticidin selection takes 7-10 days.
  5. Pick and expand at least 20 colonies. To screen the clones for those expressing the highest levels of Tet repressor, transiently transfect the positive control plasmid containing the *lacZ* gene into the cells and assay for  $\beta$ -galactosidase expression after induction with tetracycline. You will want to select for those clones exhibiting the lowest basal levels and highest inducible levels of  $\beta$ -galactosidase expression. For more information about the positive control plasmid and how to assay for  $\beta$ -galactosidase expression, please refer to the manual for the inducible expression vector that you have obtained. For more information about induction of gene expression with tetracycline, please refer to the T-REx™ System manual.
  6. Once you have obtained cell lines that stably express the Tet repressor from pcDNA6/TR, you may use these cell lines to assay for tetracycline-regulated expression of your gene of interest from the pcDNA4/TO-based expression vector.
-

# Appendix


## Blasticidin

### Blasticidin

Blasticidin S HCl is a nucleoside antibiotic isolated from *Streptomyces griseochromogenes* which inhibits protein synthesis in both prokaryotic and eukaryotic cells (Takeuchi *et al.*, 1958; Yamaguchi *et al.*, 1965). Resistance is conferred by expression of either one of two blasticidin S deaminase genes: *bsd* from *Aspergillus terreus* (Kimura *et al.*, 1994) or *bsr* from *Bacillus cereus* (Izumi *et al.*, 1991). These deaminases convert blasticidin S to a non-toxic deaminohydroxy derivative (Izumi *et al.*, 1991).

### Molecular Weight, Formula, and Structure

The formula for blasticidin S is  $C_{17}H_{26}N_8O_5 \cdot HCl$ , and the molecular weight is 458.9. The diagram below shows the structure of blasticidin.



### Handling Blasticidin

Always wear gloves, mask, goggles, and protective clothing (e.g. a laboratory coat) when handling blasticidin. Weigh out blasticidin and prepare solutions in a hood.

### Preparing and Storing Stock Solutions


Blasticidin is available from Invitrogen (Catalog no. R210-01) in 50 mg aliquots. Blasticidin is soluble in water. Sterile water is generally used to prepare stock solutions of 5 to 10 mg/ml.

- Dissolve blasticidin in sterile water and filter-sterilize the solution.
- Aliquot in small volumes suitable for one time use (see next to last point below) and freeze at  $-20^{\circ}C$  for long-term storage or store at  $+4^{\circ}C$  for short-term storage.
- Aqueous stock solutions are stable for 1-2 weeks at  $+4^{\circ}C$  and 6-8 weeks at  $-20^{\circ}C$ .
- pH of the aqueous solution should be 7.0 to prevent inactivation of blasticidin.
- Do not subject stock solutions to freeze/thaw cycles (**do not store in a frost-free freezer**).
- Upon thawing, use what you need and store the stock solution at  $+4^{\circ}C$  for up to 2 weeks.
- Medium containing blasticidin may be stored at  $+4^{\circ}C$  for up to 2 weeks.

# pcDNA6/TR Vector

## Map of pcDNA6/TR

The figure below summarizes the features of the pcDNA6/TR vector. The complete sequence for pcDNA6/TR is available for downloading from our World Wide Web site ([www.invitrogen.com](http://www.invitrogen.com)) or from Technical Service (see page 9). Please see the next page for a description of the features of the vector.



### Comments for pcDNA6/TR® 6662 nucleotides

CMV promoter: bases 232-819

Rabbit  $\beta$ -globin intron II (IVS): bases 1028-1600

*TetR* gene: bases 1684-2340

SV40 early polyadenylation sequence: bases 2346-2477

f1 origin: bases 2897-3325

SV40 promoter and origin: bases 3335-3675

EM-7 promoter: bases 3715-3781

Blasticidin resistance gene: bases 3782-4180

SV40 early polyadenylation sequence: bases 4338-4468

pUC origin: bases 4851-5521

*bla* promoter: bases 6521-6625 (complementary strand)

Ampicillin (*bla*) resistance gene: bases 5666-6526 (complementary strand)

*continued on next page*

## pcDNA6/TR Vector, continued

### Features of pcDNA6/TR

The table below describes the relevant features of pcDNA6/TR. The vector includes the rabbit  $\beta$ -globin intron II to enhance expression of the *TetR* gene. All features have been functionally tested.

Feature	Benefit
Human cytomegalovirus (CMV) immediate early promoter	Permits high-level expression of the <i>TetR</i> gene (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Nelson <i>et al.</i> , 1987)
Rabbit $\beta$ -globin intron II (IVS)	Enhances expression of the <i>TetR</i> gene (van Ooyen <i>et al.</i> , 1979)
<i>TetR</i> gene	Encodes the Tet repressor that binds to <i>tet</i> operator sequences to repress transcription of the gene of interest in the absence of tetracycline (Postle <i>et al.</i> , 1984; Yao <i>et al.</i> , 1998)
SV40 early polyadenylation signal	Permits efficient transcription termination and polyadenylation of mRNA
f1 origin	Allows rescue of single-stranded DNA
SV40 early promoter and origin	Allows efficient, high-level expression of the blasticidin resistance gene in mammalian cells and episomal replication in cells expressing SV40 large T antigen
EM-7 promoter	Synthetic prokaryotic promoter for expression of the blasticidin resistance gene in <i>E. coli</i>
Blasticidin ( <i>bsd</i> ) resistance gene	Allows selection of stable transfectants in mammalian cells (Kimura <i>et al.</i> , 1994) and transformants in <i>E. coli</i>
SV40 early polyadenylation signal	Allows efficient transcription termination and polyadenylation of mRNA
pUC origin	Permits high-copy number replication and growth in <i>E. coli</i>
<i>bla</i> promoter	Allows expression of the ampicillin ( <i>bla</i> ) resistance gene
Ampicillin ( <i>bla</i> ) resistance gene ( $\beta$ -lactamase)	Allows selection of transformants in <i>E. coli</i>

# Technical Service

---

## World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

**<http://www.invitrogen.com>**

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

---

## Contact Us

For more information or technical assistance, please call, write, fax, or email. Additional international offices are listed on our Web page ([www.invitrogen.com](http://www.invitrogen.com)).

### United States Headquarters:

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008 USA  
Tel: 1 760 603 7200  
Tel (Toll Free): 1 800 955 6288  
Fax: 1 760 602 6500  
E-mail:  
[tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)

### Japanese Headquarters

Invitrogen Japan K.K.  
Nihonbashi Hama-Cho Park Bldg. 4F  
2-35-4, Hama-Cho, Nihonbashi  
Tel: 81 3 3663 7972  
Fax: 81 3 3663 8242  
E-mail: [jpinfo@invitrogen.com](mailto:jpinfo@invitrogen.com)

### European Headquarters:

Invitrogen Ltd  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Tel (Free Phone Orders): 0800 269 210  
Tel (General Enquiries): 0800 5345 5345  
Fax: +44 (0) 141 814 6287  
E-mail: [eurotech@invitrogen.com](mailto:eurotech@invitrogen.com)

---

## MSDS Requests

To request an MSDS, please visit our Web site ([www.invitrogen.com](http://www.invitrogen.com)) and follow the instructions below.

1. On the home page, go to the left-hand column under 'Technical Resources' and select 'MSDS Requests'.
  2. Follow instructions on the page and fill out all the required fields.
  3. To request additional MSDSs, click the 'Add Another' button.
  4. All requests will be faxed unless another method is selected.
  5. When you are finished entering information, click the 'Submit' button. Your MSDS will be sent within 24 hours.
- 

*continued on next page*

## Technical Service, continued

---

### Emergency Information

In the event of an emergency, customers of Invitrogen can call the 3E Company, 24 hours a day, 7 days a week for disposal or spill information. The 3E Company can also connect the customer with poison control or with the University of California at San Diego Medical Center doctors.

3E Company  
Voice: 1-760-602-8700

---

### Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

**Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

---

# Purchaser Notification

---

## Introduction

The pcDNA6/TR vector is designed for use with the T-REX™ System available from Invitrogen. Use of the T-REX™ System and its components (“System”) is covered under a number of different licenses including those detailed below.

---

## The Nature of the Invitrogen License

Invitrogen has a license to sell the System to scientists for **academic research or one year commercial evaluation only**, under the terms described below. Use of the System for any Commercial Purpose (as defined below) other than evaluation requires the user to obtain a commercial license as detailed below. Note that such a license would cover only one part of the System. Before using the System, please read the terms and conditions set forth below. Your use of the System shall constitute acknowledgment and acceptance of these terms and conditions. If you do not wish to use the System pursuant to these terms and conditions, please contact Invitrogen’s Technical Services to return the unused and unopened System for a full credit. Otherwise, please complete the User Registration Card and return it to Invitrogen.

---

## Terms and Conditions

Invitrogen grants you a non-exclusive license to use the enclosed System for academic research or for commercial evaluation purposes only. The System is being transferred to you in furtherance of, and reliance on, such license. You may not use the System, or the materials contained therein, for any Commercial Purpose without a license for such purpose from Invitrogen.

---

## One Year Evaluation

If you are a commercial entity, your right to use the System expires after one year. Any commercial entity that wishes to use the System beyond this one-year period must obtain a commercial license from Invitrogen. Note that such a license would cover only one part of the System. Additional licenses for commercial use, as described on pages 12 and 13, may be required. Commercial entities will be contacted by Invitrogen during this one-year period regarding their desire to obtain a commercial license.

You may terminate your use of the System at any time by destroying all System components in your control. Your right to use the System will also terminate automatically if you fail to comply with the terms and conditions set forth herein. You shall, upon such termination of your rights, destroy all System components in your control, and notify Invitrogen of such in writing.

---

## Definition of Commercial Purpose

Commercial Purpose includes:

- any use of Expression Products in a Commercial Product;
- any use of Expression Products in the manufacture of a Commercial Product;
- any sale of Expression Products;
- any use (other than evaluation) of Expression Products or the System to facilitate or advance research or development of a Commercial Product; and
- any use (other than evaluation) of Expression Products or the System to facilitate or advance any research or development program the results of which will be applied to the development of Commercial Products.

“Expression Products” means products expressed with the System, or with the use of any vectors or host strains in the System. “Commercial Product” means any product intended for sale or commercial use.

---

*continued on next page*

## Purchaser Notification, continued

---

### Individual Responsibilities

Access to the System must be limited solely to those officers, employees and students of your entity who need access to perform the aforementioned research or evaluation. Each such officer, employee and student must be informed of these terms and conditions and agree, in writing, to be bound by same. You may not distribute the System or the vectors or host strains contained in it to others. You may not transfer modified, altered, or original material from the System to a third party without written notification to, and written approval from Invitrogen. You may not assign, sub-license, rent, lease or otherwise transfer any of the rights or obligations set forth herein, except as expressly permitted by Invitrogen.

---

### Governing Law

This license shall be governed in its interpretation and enforcement by the laws of the United States and California.

---

### T-REx™ System License

This product is licensed under patents pending in the U.S. and other countries and is **for research purposes only**. Inquiries about licensing for commercial or other uses should be directed to:

Licensing Coordinator  
Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008  
Phone: 760-603-7200  
Fax: 760-603-7201

---

### Use of Prokaryotic Control Elements in Eukaryotes

This product is also licensed under U.S. Patent No. 4,833,080 and corresponding patents in other countries **for research purposes only**. It may not be used for gene expression in plants nor for certain molecular modifications and screening. Inquiries about licensing for commercial manufacture or research use outside the research kits and reagents market should be directed to:

Office of Technology and Trademark Licensing  
Harvard University  
University Place, Suite 410 South  
124 Mt. Auburn Street  
Cambridge, MA 02138  
Phone: 617-495-3067  
Fax: 617-495-9568

Inquiries regarding commercial use in the research kit and reagent market should be directed to the Licensing Coordinator, Invitrogen Corporation.

---

### Blasticidin and the Blasticidin Selection Marker

Blasticidin and the blasticidin resistance gene (*bsd*) are sold under patent license and may be used for **research purposes only**. Inquiries for commercial use should be directed to:

Kaken Pharmaceutical Company, Ltd. S  
Bunkyo Green Court  
Center Office Building, 19-20 Fl.  
28-8 Honkomagome 2-chome  
Bunkyo-ku, Tokyo 113-8650, Japan  
Tel: 81 3-5977-5008  
Fax: 81 3-5977-5008

---

*continued on next page*



## Purchaser Notification, continued

---

### **Zeocin™**

Zeocin™ is a trademark of CAYLA, Toulouse, France. For commercial license information, please contact:

Licensing Coordinator  
Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008  
Phone: 760-603-7200  
Fax: 760-602-6500

---

### **CMV Promoter**

Use of the CMV promoter is covered under U.S. Patent Nos. 5,168,062 and 5,385,839 owned and licensed by the University of Iowa Research Foundation and may be used **for research purposes only**. Commercial users must obtain a license to these patents directly from the University of Iowa Research Foundation. Inquiries for commercial use should be directed to:

University of Iowa Research Foundation  
214 Technology Innovation Center  
Iowa City, IA 52242  
Phone: 319-335-4546

or

Mrs. Usha Baladrishnan  
University of Iowa Research Foundation Consultant  
Phone: 678-475-1272

---

### **BGH Polyadenylation Signal**

The BGH polyadenylation sequence is licensed under U.S. Patent No. 5,122,458 **for research purposes only**. Inquiries for commercial use should be directed to:

Bennett Cohen, Ph.D.  
Research Corporation Technologies  
101 North Wilmot Road, Suite 600  
Tucson, AZ 85711-3335  
Phone: 520-748-4400  
Fax: 520-748-0025

---

### **Product User Registration Card**

Please complete and return the enclosed Product User Registration Card for each T-REx™ System component that you purchase. This will serve as a record of your purchase and will allow Invitrogen to provide you with technical support and manual updates. It will also allow Invitrogen to update you on future developments and improvements to the T-REx™ System. The agreement outlined above becomes effective upon our receipt of your Product User Registration Card or 10 days following the sale of the T-REx™ System component to you. Use of the T-REx™ System component at any time results in immediate obligation to the terms and conditions stated in this license agreement.

---

## References

---

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989). Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. *J. Biol. Chem.* *264*, 8222-8229.
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). *Current Protocols in Molecular Biology* (New York: Greene Publishing Associates and Wiley-Interscience).
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985). A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. *Cell* *41*, 521-530.
- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Mol. Cell. Biol.* *7*, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nuc. Acids Res.* *15*, 1311-1326.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* *32*, 115-121.
- Felgner, P. L., and Ringold, G. M. (1989). Cationic Liposome-Mediated Transfection. *Nature* *337*, 387-388.
- Hillen, W., and Berens, C. (1994). Mechanisms Underlying Expression of Tn10 Encoded Tetracycline Resistance. *Annu. Rev. Microbiol.* *48*, 345-369.
- Hillen, W., Gatz, C., Altschmied, L., Schollmeier, K., and Meier, I. (1983). Control of Expression of the Tn10-encoded Tetracycline Resistance Genes: Equilibrium and Kinetic Investigations of the Regulatory Reactions. *J. Mol. Biol.* *169*, 707-721.
- Izumi, M., Miyazawa, H., Kamakura, T., Yamaguchi, I., Endo, T., and Hanaoka, F. (1991). Blasticidin S-Resistance Gene (*bsr*): A Novel Selectable Marker for Mammalian Cells. *Exp. Cell Res.* *197*, 229-233.
- Kimura, M., Takatsuki, A., and Yamaguchi, I. (1994). Blasticidin S Deaminase Gene from *Aspergillus terreus* (*BSD*): A New Drug Resistance Gene for Transfection of Mammalian Cells. *Biochim. Biophys. Acta* *1219*, 653-659.
- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987). Negative and Positive Regulation by a Short Segment in the 5'-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. *Mol. Cell. Biol.* *7*, 4125-4129.
- Postle, K., Nguyen, T. T., and Bertrand, K. P. (1984). Nucleotide Sequence of the Repressor Gene of the Tn10 Tetracycline Resistance Determinant. *Nuc. Acids Res.* *12*, 4849-4863.
- Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* *6*, 742-751.
- Takeuchi, S., Hirayama, K., Ueda, K., Sakai, H., and Yonehara, H. (1958). Blasticidin S, A New Antibiotic. *The Journal of Antibiotics, Series A* *11*, 1-5.
- van Ooyen, A., van den Berg, J., Mantei, N., and Weissmann, C. (1979). Comparison of Total Sequence of a Cloned Rabbit Beta-globin gene and its Flanking Regions With a Homologous Mouse Sequence. *Science* *206*, 337-344.

---

*continued on next page*

## References, continued

---

Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. *Cell* *11*, 223-232.

Yamaguchi, H., Yamamoto, C., and Tanaka, N. (1965). Inhibition of Protein Synthesis by Blastocidin S. I. Studies with Cell-free Systems from Bacterial and Mammalian Cells. *J. Biochem. (Tokyo)* *57*, 667-677.

Yao, F., Svensjö, T., Winkler, T., Lu, M., Eriksson, C., and Eriksson, E. (1998). Tetracycline Repressor, tetR, Rather than the tetR-Mammalian Cell Transcription Factor Fusion Derivatives, Regulates Inducible Gene Expression in Mammalian Cells. *Hum. Gene Ther.* *9*, 1939-1950.

---

©1999-2002 Invitrogen Corporation. All rights reserved.